

Section II: Preparation of Agarose Gels

Suggested Agarose Concentrations and Dye Migration Information

Table 1: Suggested agaroses for particular applications

Size Range (Base Pairs)	Agarose Type	Application
20 – 800	MetaPhor® Agarose	High resolution analysis; 2% size differences
50 – 1,000	NuSieve® 3:1 Agarose	Analysis and blotting; 4% – 6% size differences resolved
	NuSieve® GTG® Agarose	Analysis and blotting; In-gel; 6% size differences resolved
1,000 – 10,000	SeaKem® GTG® Agarose	Analysis and blotting; recovery required
	SeaPlaque® GTG® Agarose	In-gel
10,000	SeaKem® Gold Agarose	Analysis

Table 3: Suggested agarose concentrations for DNA sizes

Size Range (Base Pairs)	Final Agarose Concentration % (w/v)	
	1X TAE Buffer	1X TBE Buffer
SeaKem® LE and SeaKem® GTG® Agarose		
1,000 – 23,000	0.60	0.50
800 – 10,000	0.80	0.70
400 – 8,000	1.00	0.85
300 – 7,000	1.20	1.00
200 – 4,000	1.50	1.25
100 – 3,000	2.00	1.75
NuSieve® 3:1 Agarose		
500 – 1,000	3.0	2.0
100 – 500	4.0	3.0
10 – 100	6.0	5.0
MetaPhor® Agarose		
150 – 800	2.0	1.8
100 – 600	3.0	2.0
50 – 250	4.0	3.0
20 – 130	5.0	4.0
<80	—	5.0
SeaPlaque® and SeaPlaque® GTG® Agarose		
500 – 25,000	0.75	0.70
300 – 20,000	1.00	0.85
200 – 12,000	1.25	1.00
150 – 6,000	1.50	1.25
100 – 3,000	1.75	1.50
50 – 2,000	2.00	1.75
NuSieve® GTG® Agarose		
500 – 1,000	2.5	2.0
150 – 700	3.0	2.5
100 – 450	3.5	3.0
70 – 300	4.0	3.5
10 – 100	4.5	4.0
8 – 50	5.0	4.5
SeaKem® Gold Agarose†		
5,000 – 50,000	0.3	—
1,000 – 20,000	0.5	—
800 – 10,000	0.8	—
400 – 8,000	1.0	—

† TBE Buffer is not recommended for separation of DNA >12,000 bp.

Table 2: Properties of TAE and TBE Buffer Systems

Buffer	Suggested Uses and Comments
TAE Buffer	Use when DNA is to be recovered Use for electrophoresis of large (>12 kb) DNA Low ionic strength Low buffering capacity – recirculation may be necessary for extended electrophoretic times
TBE Buffer	Use for electrophoresis of small (<1 kb) DNA Decreased DNA mobility High ionic strength High buffering capacity – no recirculation required for extended run times

Table 4: Migration of double-stranded DNA in relation to Bromophenol Blue (BPB) and Xylene Cyanol (XC) in agarose gels

1X TAE Buffer			1X TBE Buffer		
XC	BPB	% Agarose	XC	BPB	
SeaKem® LE and SeaKem® GTG® Agarose					
24,800	2,900	0.30	19,400		2,850
11,000	1,650	0.50	12,000		1,350
10,200	1,000	0.75	9,200		720
6,100	500	1.00	4,100		400
3,560	370	1.25	2,500		260
2,800	300	1.50	1,800		200
1,800	200	1.75	1,100		110
1,300	150	2.00	850		70
NuSieve® 3:1 Agarose					
950	130	2.50	700		70
650	80	3.00	500		40
350	40	4.00	250		20
200	30	5.00	140		8
120	20	6.00	90		4
MetaPhor® Agarose					
480	70	2.00	310		40
200	40	3.00	140		35
120	35	4.00	85		30
85	30	5.00	60		15
SeaPlaque® and SeaPlaque® GTG® Agarose					
11,700	1,020	0.50	6,100		400
4,000	500	0.75	2,850		280
2,300	350	1.00	1,700		180
1,500	200	1.25	1,000		100
1,000	150	1.50	700		70
700	100	1.75	500		50
550	60	2.00	400		30
320	30	2.50	250		10
NuSieve® GTG® Agarose					
750	175	2.50	460		75
400	120	3.00	210		35
115	<20	4.00	150		<20
100	<20	5.00	80		<20
85	<20	6.00	50		<20
SeaKem® Gold Agarose					
24,800	3,550	0.30	19,000		2,550
12,200	2,050	0.50	9,200		1,500
9,200	1,050	0.75	7,100		800
6,100	760	1.00	4,000		500
4,100	600	1.25	2,550		350
2,600	400	1.50	1,900		250
2,000	330	1.75	1,400		180
1,500	250	2.00	1,000		100