BD Vacutainer®

CPT™ Mononuclear Cell Preparation Tube

Helping all people
live healthy lives

One-step cell separation

BD Vacutainer®

CPT™ Mononuclear Cell Preparation Tube

The *only* one-step, closed system for collection, mononuclear cell separation, and transportation


The BD Vacutainer® CPT™ tube is a one-step, closed system for blood collection, mononuclear cell separation, and transportation. It is an evacuated tube containing an anticoagulant (either sodium citrate or sodium heparin) and a cell separation medium comprised of a polyester gel and a density gradient liquid. The polyester gel forms a stable barrier between the cell layers after centrifugation to allow for transportation to the site where testing is to be performed.

All this – achieved in one tube!


Before Blood Collection

After Centrifugation

Mononuclear Cell Separation

One-step Cell Separation

- › Separate mononuclear cells in a single step from whole blood
- › Standardize cell separation protocol
- › Obtain high mononuclear cell viability and yield
- › Transport separated cells in the same processed tube
- › Reduce risk of specimen contamination
- › Reduce biohazard exposure
- › For *In Vitro* Diagnostic use

It's easy to Process BD Vacutainer® CPT™ Tubes

1

Blood Collection

Use standard venipuncture procedure.

2

Centrifugation

Gently invert 8 times before centrifugation at room temperature:
Heparin: 15 min.
Citrate: 20 min.
Horizontal rotor with swing-out head and appropriate tube adapters at 1500xg RCF.


Note: For best results, centrifuge within 2 hours after blood collection, but no later than 24 hours after collection!

Acceptable blood volume:
8 mL
7 mL
6 mL

Anticoagulated Whole Blood

Gel Barrier

Density Gradient Fluid


3

Result

Mononuclear cell separation is completed in one step, in the same tube!

Plasma
Lymphocytes & Monocytes
Gel Barrier
Erythrocytes & Neutrophils

Density Gradient Fluid


4

Transportation

To re-suspend cells, gently invert one time before transportation.

Cell Suspension (plasma and mononuclear cells)

BD Vacutainer®

CPT™ Mononuclear Cell Preparation Tube

Product Description

The BD Vacutainer® CPT™ tube is a closed, one-step system for blood collection, mononuclear cell separation and transportation. It is an evacuated tube containing an anticoagulant (either Sodium Citrate or Sodium Heparin) and a cell separation medium comprised of a polyester gel and a density gradient liquid.

The polyester gel forms a stable barrier between the cell layers after centrifugation to allow for transportation to the site where the assay is to be performed.

Unique Benefits

- Separate mononuclear cells in a single step from whole blood
- Obtain high cell viability and yield
- Transport separated cells in the same processed tube
- Reduce risk of specimen contamination and biohazard exposure
- For *In Vitro* Diagnostic use

« BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company »

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