
With a variety of options to choose from, **SBI accelerates exosome isolation** with products that enable fast, multiplexed exosome enrichment from a wide range of biofluids.

SALIVA<sup>1</sup> / URINE<sup>2</sup> / FOLLICULAR FLUID<sup>3</sup> / PLASMA<sup>4</sup> / SERUM<sup>5</sup> / TISSUE CULTURE MEDIA<sup>6</sup> / BREAST MILK<sup>7</sup> / ASCITES FLUID<sup>8</sup> / MORE

## ExoQuick Family

The ExoQuick portfolio of products originated with the ExoQuick<sup>®</sup> and ExoQuick-TC<sup>®</sup> isolation kits that enable high-throughput, high-yield exosome isolation. Using these unique, polymer-based reagents, exosomal vesicles are gently and reliably precipitated from solution. Compatible with virtually any biofluid (such as serum, plasma, tissue culture media, etc) and a wide variety of downstream applications, ExoQuick and ExoQuick-TC are effective and proven<sup>3,4,9</sup> alternatives to ultracentrifugation.

Quick and easy—simple protocol isolates exosomes in as little as 30 minutes

Universal—enables exosome isolation from all biofluids tested to-date

Broadly usable—isolated exosomes perform well in a variety of downstream applications, including miRNA profiling, NGS, and mass spec analysis

Functional—isolated exosomes are intact and bioactive for functional studies and engineering



Exosomal and Non-Exosomal Transport of Extra-Cellular microRNAs in Follicular Fluid: Implications for Bovine Oocyte Developmental Competence. Sohel, M. M. H. et al. PLoS ONE 8, (2013). PMID: PMC3817212

## ExoQuick ULTRA

Drawing upon years of exosome experience, the SBI team has pushed ExoQuick extracellular vesicle isolation technology to new peaks of performance with ExoQuick ULTRA for Serum and Plasma and ExoQuick-TC ULTRA for Tissue Culture Media and other biofluids. While many EV isolation methods require a choice between high yields, easy protocols, clean preps, and low costs, ExoQuick ULTRA is able to deliver on all fronts for trade-off free EV preparation.

- Even cleaner—reduces carry-over of albumins and immunoglobulins compared to ultracentrifugation
- Higher yields—isolate more EVs per normalized input volume than UC and other kits
- Better biomarker detection—see what you've been missing when you increase the sensitivity of EV biomarker detection
- Fast—requires < 20-minutes of hands-on time
- Cost-effective—save money with each reaction compared to using competitor kits

⚡ Fluorescent nanoparticle tracking analysis demonstrates the high EV yields delivered by ExoQuick ULTRA compared to UC. Comparison of different isolation methods on EV yields by both volume of input serum (per mL, top panel) and amount of isolated protein (per mg as measured by fluorometric Qubit protein assay, bottom panel).

⚡ ExoQuick ULTRA delivers high yields of clean exosomes from serum. Western blot shows that the ExoQuick ULTRA prep contains the highest levels of exosome-specific markers CD9, CD81, and Hsp70 and the lowest level of total carryover protein.

### 🔍 VALIDATION DATA

Find additional validation data and references using ExoQuick ULTRA—visit [systembio.com/exoquick-ultra](https://systembio.com/exoquick-ultra)

🔬 Biochemical and biologic characterization of exosomes and microvesicles as facilitators of HIV-1 infection in macrophages. Kadiu, I. et al. J Immunol. 2012 Jul 15; 189(2):744-54. PMID: PMC3786185

### 💡 CONSIDER THIS

With just 100 µl plasma, ascites fluid, or serum, you can isolate exosomes from small animal models for biomarker analysis.

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