
Small RNA sequencing library preparation from total RNA and biofluids


RealSeq[®]-AC: Reducing bias in small-RNA sequencing

- Accurately quantifies biologically relevant small RNAs
- Eliminates bias-induced over/under represented small RNAs
- Allows discovery of novel small RNAs
- Inputs between 1 ng to 1 µg of total RNA


One pmole of miRXPlore Universal Reference pool (Miltenyi Biotec) was used to compare incorporation bias in six different commercially available library preparation kits. Purified libraries were sequenced on the Illumina MiSeq platform. Trimmed sequencing reads were aligned to a custom miRNA reference. Reads mapping to miRNAs were counted and fold-deviations from the equimolar input were calculated and plotted as log2 values. Measurements of miRNA levels within a factor of two of the expected values (between vertical lines) are considered unbiased (Fuchs et al, 2015). The method of adapter attachment to the miRNAs is noted in the legend.

RealSeq[®]-AC sequence-ready library in one day

Small RNA sequencing library preparation from total RNA and biofluids

RealSeq[®]-AC schematic


RealSeq[®]-AC is highly efficient, detecting more miRNAs in total RNA samples

	Kit I	Kit N	Kit B	Kit T	RealSeq [®] -AC
miRNAs with >5 reads	404	452	412	324	500
miRNAs with >10 reads	328	365	352	239	385

RealSeq[®]-AC detects more miRNAs with over 5 or 10 reads per million, respectively, from total RNA samples compared to other kits for miRNA sequencing library preparation. RealSeq[®]-AC is optimized for inputs between 1 ng to 100 ng of total RNA. Lower numbers of PCR cycles are reducing PCR-induced issues.

Elimination of biases using RealSeq® - Methods comparison from Barberan-Soler et al., 2018*

Results obtained with the miRXPlore Universal Reference pool (963 miRs)


- RealSeq®-AC detects 71,8% of miRs without bias
- Other kits detect only between 12,9% and 48% of miRs

* Decreasing miRNA sequencing bias using a single adapter and circularization approach. Sergio Barberán-Soler, Jenny M. Vo, Ryan E. Hogans, Anne Dallas, Brian H. Johnston and Sergei A. Kazakov. Genome Biology 2018 19:105


[Lire la publication](#)

Detection of a large diversity of miRs in biofluids

Percentage reads for each plasma miRNA with Kit Q

Percentage reads for each plasma miRNA with RealSeq-Biofluids

Due to low sequencing biases, RealSeq®-biofluids allow more accurate analysis of biological samples

- Kit Q identifies miR-16-5p as the most abundant plasma miRNA with 43.1% of all reads
- However; miRXPlore pool sequencing showed that Kit Q over-represents this miRNA by >5 fold
- In contrast RealSeq®-biofluids identifies miR-221-3p as the the most abundant miRNA with 11.8% of all reads
- miRXPlore pool sequencing showed that RealSeq® detects this miRNA without biases
- In both experiments, miRs purification has been carried out with Quick-cfRNA Serum & Plasma Kit (Zymo)

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	Référence	Désignation	Conditionnement
A partir d'ARN total 1 ng à 1 µg Billes de purification non incluses*	REA500-00012*	RealSeq-AC small RNA Library kit for Illumina sequencing	12 libraries
	REA500-00024*	RealSeq-AC small RNA Library kit for Illumina sequencing	24 libraries
	REA500-00048*	RealSeq-AC small RNA Library kit for Illumina sequencing	48 libraries
A partir de cfRNA isolés de 50 µl de plasma ou sérum Billes de purification incluses	REA600-00012	RealSeq-Biofluids small RNA Library kit for Illumina sequencing	12 libraries
	REA600-00024	RealSeq-Biofluids small RNA Library kit for Illumina sequencing	24 libraries
	REA600-00048	RealSeq-Biofluids small RNA Library kit for Illumina sequencing	48 libraries

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Billes de purification Beckman	REA510-00012	RealSeq-AC size select beads	12 libraries
	REA510-00024	RealSeq-AC size select beads	24 libraries
	REA510-00048	RealSeq-AC size select beads	48 libraries

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