


# *in vivo* transfection

## *in vivo*-jetPEI®

- ⊕ Polyvalent: *in vivo* delivery of DNA, si/sh/miRNA in any animal model
- ⊕ Easy to use: Two-step protocol
- ⊕ Renowned: Most advanced *in vivo* transfection technology with over 700 peer-reviewed publications
- ⊕ Successful: Used from fundamental research to Human clinical trials




### *in vivo*-jetPEI® can be used in a wide range of organs:



Complexes were formed following standard protocols.

### Easy-to-use 2 steps protocol:



**Nous contacter**

| Product                 | Part N°   | Reagent size | Amount of 10% glucose solution |
|-------------------------|-----------|--------------|--------------------------------|
| <i>in vivo</i> -jetPEI® | 101000040 | 0.1 mL       | 10 mL                          |
|                         | 101000030 | 0.5 mL       | 2 x 10 mL                      |