


# Cellometer® Auto 1000

## Cell Counter

### Accurate Cellometer Cell Counting ... Now Even Faster and Easier

- ➔ Improve counting of live/dead cells, individual cells within clusters, and irregular-shaped cells: Proprietary Cellometer software
- ➔ Spend less time counting and more time on research: Concentration and viability in <10 seconds!
- ➔ Easily count cells in the hood: Stand-alone system with a large, user-friendly touch screen
- ➔ View, print, and save counted cell images, data reports, and cell size histograms: Connect to a network or use a USB drive




## Simple, Automated Cell Counting

Improve Speed, Accuracy, and Consistency



### Load 20 µl of Sample




### Insert Counting Chamber



### Select Assay & Click Count



### View Image & Results



Tutoriel Utilisation Cellometer®  
Auto 1000

**Service Technique** : Réactifs : 01 34 60 60 24 – tech@ozyne.fr - Instrumentation : 01 30 85 92 89 - equipement@ozyne.fr  
**Service Commande - Clients** : 01 34 60 15 16 – commande@ozyne.fr  
Nous contacter : [www.ozyne.fr/contact](http://www.ozyne.fr/contact)