

MycoAlert™ PLUS Assay

The 20 Minute Route to Mycoplasma Detection



MycoAlert™ PLUS Assay
for Mycoplasma Detection



For a decade, Lonza's MycoAlert™ Assay has been known as an easy-to-use bioluminescence-based mycoplasma test, providing reliable results in 20 minutes. It exploits the activity of enzymes found in the vast majority of mycoplasma species.


MycoAlert™ PLUS Assay

This improved assay has a 20 times higher light output, thus expanding compatibility of the established technology to wider range of instruments and samples.

Contact your local sales representative for more information.

Benefits:

- No DNA extraction required, thanks to bioluminescence-based technology
- Simple process with just two pipetting steps
- MycoAlert™ PLUS Assay for use with less sensitive plate luminometers or multifunctional readers
- MycoAlert™ PLUS Assay for testing of unused media, media supplements or water



Performance of MycoAlert™ PLUS Assay compared to MycoAlert™ Assay. A dilution series of MycoAlert™ Assay Control demonstrates the increased sensitivity of various luminometer models when using MycoAlert™ PLUS Assay, compared to the first generation MycoAlert™ Assay.

Ordering Information

Cat. No.	Description	Quantity
LT07-701	MycoAlert™ PLUS Mycoplasma Detection Kit	10 tests
LT07-703	MycoAlert™ PLUS Mycoplasma Detection Kit	30 tests
LT07-705	MycoAlert™ PLUS Mycoplasma Detection Kit	50 tests
LT07-710	MycoAlert™ PLUS Mycoplasma Detection Kit	100 tests
LT07-518	MycoAlert™ Assay Control	10 tests



Lonza Walkersville, Inc. – Walkersville, MD 21793

For research use only. Not for use in diagnostic procedures. Unless otherwise noted, all trademarks herein are marks of the Lonza Group or its affiliates. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.

© Copyright 2012, Lonza Walkersville, Inc. All rights reserved.
FL-Myco20min 08/12 CD-TS044

Nous contacter



Service Technique (réactifs et instrumentation) : 01 34 60 60 24 - tech@ozyme.fr
Service Commande - Clients : 01 34 60 15 16 – commande@ozyme.fr