

COMPATIBILITÉ APPAREILS ET PREMIX DE QPCR

	Applied Biosystems 7000/7700 - 7300 - StepOne - 7900HT	Applied Biosystems 7500 - 7500 Fast - ViiA™ 7 - QuantStudio™	LightCycler® (Roche Diagnostics)	Smart Cycler® (Cepheid)	CFX96™ - iCycler® - Chromo4™ - Opticon™ (Bio-RAD)	Rotorgen (Qiagen ou Corbette)	EP RealPlex (Eppendorf)	AiraMX - Mx3000P - Mx300x (Agilent)	Eco qPCR Illumina
ONE Green® Fast qPCR Premix*	OZYA008	OZYA008	OZYA008	OZYA008	OZYA008	OZYA008	OZYA008	OZYA008	OZYA008
SensiFAST™ SYBR®§	Hi-ROX	Lo-ROX	No-ROX	No-ROX	No-ROX	No-ROX	No-ROX	Lo-ROX	No-ROX
SensiFAST™ Probe§	Hi-ROX	Lo-ROX	No-ROX	No-ROX	No-ROX	No-ROX	No-ROX	Lo-ROX	No-ROX
Quel ROX ?	ROX I (high concentration)	ROX II (low concentration)	no ROX	no ROX	no ROX	no ROX	no ROX	ROX II (low concentration)	no ROX

* ROX I et ROX II fournis dans des tubes séparés

§ ROX dye intégré dans le premix