

Comment choisir son tampon de lyse ?

Type de protéine d'intérêt Type d'échantillon	Protéine totale			Protéine phosphorylée				
	Tampons	Ajout(s) recommandé(s)		Tampons	Ajout(s) recommandé(s)			
Cellules	lysis buffer #9803	Protease Inhibitor Cocktail (100X) #5871	PMSF #8553	lysis buffer #9803	Protease Inhibitor Cocktail (100X) #5871	PMSF #8553	Phosphatase Inhibitor Cocktail (100X) #5870	Protease/Phosphatase Inhibitor Cocktail (100X) #5872
Tissus	RIPA buffer #9806	Protease Inhibitor Cocktail (100X) #5871	PMSF #8553	RIPA buffer #9806	Protease Inhibitor Cocktail (100X) #5871	PMSF #8553	Phosphatase Inhibitor Cocktail (100X) #5870	Protease/Phosphatase Inhibitor Cocktail (100X) #5872