

Neuronal and Glial Cell Marker Atlas

Cell Type Markers

*Proliferating Cells

Neural Stem Cell Markers

- Musashi
- Notch1
- SSEA1

Neuroepithelial Cells*

- E-Cadherin
- HES1
- Nestin
- Notch1
- Sox2
- Sox10

Schwann Cells

- EAAT1 (GLAST-1)
- EAAT2 (GLT-1)
- GAP43
- GFAP
- Myelin Basic Protein
- NCAM
- p75NTR
- S100B

Radial Glia*

- BLP (FABP7)
- N-Cadherin
- GFAP
- EAAT1 (GLAST-1)
- Nestin
- Sox2
- Tenascin C
- Vimentin

Microglia



- ASC
- CD11b (ITGAM)
- CD40
- CD45
- CD68
- F4/80
- Iba1
- PU.1
- TMEM119

Homeostatic:

- C1qA/B
- CX3CR1
- Cystatin 3 (CST3)
- P2RY12

Functional Markers:

- Arg3.1 (Arc)
- CaMKII
- CREB
- EGR1
- Erk
- c-Fos
- PKA

Proliferating Cell Markers:

- BrdU
- Ki-67
- PCNA

Oligodendrocytes

- A2B5
- CNPase
- MAG
- MOG
- Myelin Basic Protein
- NG2
- Olig2
- SOX10

Astrocytes

- A2B5
- ALDH1L1
- EAAT1 (GLAST-1)
- EAAT2 (GLT-1)
- GFAP
- Notch
- S100B

Intermediate Progenitors*

- EOMES (Tbr2)
- Neurogenin 2
- Pax6

Immature Neurons

- Doublecortin
- NeuroD1
- TBR1
- β 3-Tubulin

Mature Neurons

- MAP2
- NeuN
- Neurofilament-H
- Neurofilament-L
- Neurofilament-M
- Neuron Specific Enolase (NSE)
- Tau
- Thy1
- β 3-Tubulin
- UCHL1

Neurotransmitters

Neuronal Subtypes

Glutamatergic Neurons

- VGLUT1/2

GABAergic Neurons

- ADORA2A
- DARPP-32
- GAD1/2
- Neuropeptide Y
- Penk

Dopaminergic Neurons

- ALDH1A1
- FoxA2
- GIRK2
- Lmx1B
- Tyrosine Hydroxylase

Serotonergic Neurons

- FEV
- SERT
- TPH2
- Tryptophan Hydroxylase

Cholinergic Neurons

- ChAT

Interneurons

- Calbindin
- Calretinin
- 5-HT3A Receptor
- Parvalbumin
- Somatostatin
- VIP

Motor Neurons

- ChAT
- HB9
- Islet-1/2
- Neurogenin 2
- Olig2

Neurodegenerative Disease-associated Markers

Alzheimer's Disease

- β -Amyloid (A β)
- ApoE
- CD33
- Presenilin 1/2 (PSEN1/2)
- Tau
- TREM2


Parkinson's Disease

- DJ-1
- GBA
- LRRK2
- Parkin
- PINK1
- α -Synuclein
- UCHL1
- VPS35

ALS

- C9orf72
- FUS
- SOD1
- SQSTM1
- TDP43


Disease-associated Microglia (DAM)



Stage 1

- Upregulated
- ApoE
- DAP12 (TYROBP)
- TREM2

- Downregulated
- P2RY12
- CX3CR1
- TMEM119
- CD33



Stage 2 (TREM2-dependant)

- Axl
- Cystatin F (CST7)
- Dectin-1 (CLEC7a)
- CD11c (ITGAX)
- Spp1

Key Products for Neuronal and Glial Cell Markers


Visit www.cellsignal.com for a full listing of antibodies, related kits, and reagents.

ANTIBODIES & KITS


	APPLICATIONS	REACTIVITY
#14974 β -Amyloid (1-42 Specific) (D9A3A) Rabbit mAb	WB, IF-F	H, (M, R)
#15126 β -Amyloid (D3D2N) Mouse mAb	WB, IHC-P, IF-F	H
#8243 β -Amyloid (D54D2) XP® Rabbit mAb	WB, IP, IF-P	H
#13366 ApoE (pan) (D7I9N) Rabbit mAb	WB, IP, IHC-P	H
#13176 Calbindin (D114Q) XP® Rabbit mAb	WB, IP, IHC-P, IF-F	H, M, R
#40868 CD40 (D8W3N) Rabbit mAb	WB, IHC-P	H
#13917 CD45 (D9M8I) XP® Rabbit mAb	WB, IP, IHC-P, IF-IC, F	H
#76437 CD68 (D4B9C) XP® Rabbit mAb	IHC-P, IF-IC, F	H
#5664 CNPase (D83E10) XP® Rabbit mAb	WB, IP, IF-F	H, M, R
#9197 CREB (48H2) Rabbit mAb	WB, IP, IHC-P, IHC-F, IF-F, IF-IC, F, ChIP, ChIP-seq	H, M, R, Mk, Dm
#2306 DARPP-32 (19A3) Rabbit mAb	WB, IP, IHC-P, IF-F	M, R, (H)
#14802 Doublecortin (A8L1U) Rabbit mAb	IF-F	M
#3195 E-Cadherin (24E10) Rabbit mAb	WB, IHC-P, IHC-F, IF-IC, F	H, M, (B, Dg, Pg)
#14472 E-Cadherin (4A2) Mouse mAb	WB, IP, IHC-P, IF-IC, F	H, M, R
#5684 EAAT1 (D44E2) XP® Rabbit mAb	WB, IP, IF-F	H, M, R
#3838 EAAT2 Antibody	WB, IF-F	M, R, (H)
#70076 F4/80 (D2S9R) XP® Rabbit mAb	WB, IP, IHC-P	M
#8186 FoxA2/HNF3 β (D56D6) XP® Rabbit mAb	WB, IP, IF-IC	H, M, R
#8945 GAP43 (D9C8) Rabbit mAb	WB, IP, IF-F	H, M, R
#12389 GFAP (D1F4Q) XP® Rabbit mAb	WB, IF-F	H, M, R
#3670 GFAP (GA5) Mouse mAb	WB, IHC-P, IF-F, IF-IC	H, M, R
#49420 ITGAM/CD11B (D6X1N) Rabbit mAb	WB, IHC-P	H
#9449 Ki-67 (8D5) Mouse mAb	IHC-P, IF-IC, F	H
#9129 Ki-67 (D3B5) Rabbit mAb	IF-F, IF-IC, F	H, M, R
#9043 MAG (D4G3) XP® Rabbit mAb	WB, IP, IF-F	H, M, R
#8707 MAP2 (D5G1) XP® Rabbit mAb	WB, IP, IF-F	M, R, (H)
#78896 Myelin Basic Protein (D8X4Q) XP® Rabbit mAb	WB, IHC-P, IF-F, IF-IC	H, M, R
#14215 N-Cadherin (13A9) Mouse mAb	WB, IP, IF-IC	H, M, R, Mk
#13116 N-Cadherin (D4R1H) XP® Rabbit mAb	WB, IP, IHC-P, IF-IC	H, M
#3576 NCAM (CD56) (123C3) Mouse mAb	WB, IHC-P, IF-IC, F	H
#4760 Nestin (Rat-401) Mouse mAb	IHC-P, IF-F	R
#24307 NeuN (D4G4Q) XP® Rabbit mAb	WB, IHC-P, IF-F	H, M, R
#2836 Neurofilament-H (RMdO 20) Mouse mAb	WB, IHC-P, IF-F	H, M, R
#2837 Neurofilament-L (C28E10) Rabbit mAb	WB, IHC-P, IF-F	H, M, R
#2835 Neurofilament-L (DA2) Mouse mAb	WB, IHC-P, IF-F	H, M, R
#2838 Neurofilament-M (RMdO 14.9) Mouse mAb	WB, IP, IHC-P, IF-F	H, M, R
#11976 Neuropeptide Y (D7Y5A) XP® Rabbit mAb	WB, IHC-P, IF-F, IF-IC	H, M, R
#3608 Notch1 (D1E11) XP® Rabbit mAb	WB, IP, IHC-P, ChIP	H, M, R
#8238 p75NTR (D4B3) XP® Rabbit mAb	WB, IP, IF-IC	H, M, R
#13110 PCNA (D3H8P) XP® Rabbit mAb	WB, IP, IHC-P, IF-F, IF-IC, F	H, M, R, Mk
#36233 PSD95 (7E3) Mouse mAb	WB, IF-F, IF-IC	H, M, R
#3450 PSD95 (D27E11) XP® Rabbit mAb	WB, IF-F	H, M, R
#89356 Sox10 (D5V9L) Rabbit mAb	WB, IP, IF-IC	H

ANTIBODIES & KITS


	APPLICATIONS	REACTIVITY
#3579 Sox2 (D6D9) XP® Rabbit mAb	WB, IF-IC, F	H, (B, Dg, Hr)
#4900 Sox2 (L1D6A2) Mouse mAb	WB, IF-IC, F	H, M, (R, B, Dg, Hr)
#5297 Synapsin-1 (D12G5) XP® Rabbit mAb	WB, IP, IHC-P, IF-F	H, M, R
#4179 α -Synuclein (D37A6) XP® Rabbit mAb	WB, IP, IHC-P, IF-F	M, R
#5568 β 3-Tubulin (D71G9) XP® Rabbit mAb	WB, IP, IF-F, IF-IC	H, M, R
#4466 β 3-Tubulin (TU-20) Mouse mAb	WB, IHC-P, IF-F	H, M, R
#46687 Tau (D1M9X) XP® Rabbit mAb	WB, IHC-P, IF-F, IF-IC	H, M, R
#12885 Phospho-Tau (Thr181) (D9F4G) Rabbit mAb	WB, IP, IHC-P	H, M, R
#3448 TDP43 (G400) Antibody	WB, IP, IF-IC	H, M, R
#91068 TREM2 (D8I4C) Rabbit mAb	WB, IP, IF-IC	H
#13106 Tyrosine Hydroxylase (A8Y7R) Rabbit mAb	WB, IF-F, IF-IC	H, M, R
#13179 UCHL1 (D3T2E) XP® Rabbit mAb	WB, IHC-P, IF-IC, F	H, M, R, Mk, (Hm, B, Dg, Pg, Hr)
#5741 Vimentin (D21H3) XP® Rabbit mAb	WB, IHC-P, IF-IC, F	H, M, R, Mk




F4/80 (D2S9R) XP® Rabbit mAb
#70076: IHC analysis of paraffin-embedded mouse liver using #70076.



α -Synuclein (D37A6) XP® Rabbit mAb
#4179: IHC analysis of paraffin-embedded mouse brain using #4179.



GFAP (D1F4Q) XP® Rabbit mAb
#12389: Confocal IF analysis of rat cerebellum using #12389 (green) and Neurofilament-H (RMdO 20) Mouse mAb #2836 (red). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye).



β -Amyloid (D54D2) XP® Rabbit mAb
#8243: Confocal IF analysis of paraffin-embedded human Alzheimer's brain using #8243 (green) and Tau (Tau46) Mouse mAb #4019 (red). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye).

APPLICATIONS & REACTIVITY KEY: WB Western Blotting / IP Immunoprecipitation / IHC Immunohistochemistry / IF Immunofluorescence / F Flow Cytometry / ChIP Chromatin Immunoprecipitation / -IC Immunocytochemistry / -P Paraffin / -F Frozen / E-P Peptide ELISA / H human / M mouse / R rat / Hm hamster / Mk monkey / C chicken / Mi mink / Dm D. melanogaster / X Xenopus / Z zebrafish / B bovine / Dg dog / Sc S. cerevisiae / All all species expected / () 100% sequence homology

© 2018 Cell Signaling Technology, Inc. Cell Signaling Technology, CST, and XP are trademarks of Cell Signaling Technology, Inc. U.S. Patent No. 5,675,063. All other trademarks are the property of their respective owners.

For pathway key and background please visit:

www.cellsignal.com/neuro-atlas

